

## Automated Pensky-Martens Closed Cup Flash Point Tester approx 100

# Simple Reliable Robust



- Easy Maintenance
- Enhanced Safety Functions
- Customized Options
- User-Friendly

#### **Easy Replacement and Maintenance**

 One-touch arm cover removal for easier igniter replacement and maintenance.



#### **Safety Functions**

- Overcut (overheat prevention)
- Fault diagnosis
- ●CO<sub>2</sub> / N<sub>2</sub> fire containment (external connection option)

#### **Emergency Switch**

Additional emergency switch at front panel for increased safety.

#### **Customized Options**

- 7-inch color LCD touch screen
- Flexible Custom modes to allow free setting of stirring speed, heating rate, ignition source application timing, pre-heating of sample, etc.

#### I/O Check Test operation and fault diagnosis can be performed. Ready Menu > System > 1/0 Check Photo Sensor (Application) ThermoFuse Connection CLOSE OPEN OFF ON Photo Sensor (Arm Position) Fire Detection OFF ON OFF ON 1000 Application Ignition 100 Stirrer 2023/10/30 1325

**Specifications** 

C€ KK


Product Name	Automated Pensky-Martens Closed Cup Flash Point Tester apm-100
Test Methods	ASTM D93, ISO 2719, ISO 15267, IP 34, JIS K 2265-3
Measuring Range	Ambient to 370 °C
Test Modes	Preset Test Modes: ASTM, JIS, and ISO (Compliant with Various Test Procedures) Custom Test Modes: Customizable Modes for Various Testing Conditions
Display	7 inch touchscreen color LCD display
Temperature Sensor	PT-100 in stainless steel sheath
Flash Detection	Temperature Detection (CRC thermocouple)
Ignition Source	Gas ignition with automatic lighting or Electric ignition. Interchangeable.
Barometric Pressure Correction	Barometric pressure correction automatic calculation
Safety Functions	<ul> <li>①Overcut (overheat prevention)</li> <li>Expected flash point +20 °C (No flash)</li> <li>370 °C (No flash)</li> <li>②Fault diagnosis</li> <li>Disconnection of Temperature Sensor</li> <li>Disconnection of Flash Detector</li> <li>Disconnection of Igniter etc.</li> <li>③CO<sub>2</sub> / N<sub>2</sub> fire containment (external connection option)</li> </ul>
External Output	USB x 2, Ethernet x 1
Power Supply	AC100 V to 120 V / AC220 V to 240 V 50 / 60 Hz 900 VA
Operating Environment	Temperature: 10 °C to 35 °C Humidity: 80 % RH or less (no condensation)
Dimensions	W 240 mm x D 480 mm x H 400 mm
Weight	18 kg
Minimum Footprint	W 340 mm x D 580 mm x H 700 mm


TANAKA reserves the right of changes without prior notice.



Tanaka File Export System (TFES) is a PC software that collects test data from Tanaka instruments and saves it to a designated folder as a text or CSV file.

If a LIMS or other data management systems are set up to collect the data from the folder, the data can be incorporated into LIMS.





#### **Standard Accessories**

Code No.	Part Name	Qty	Remarks
1007478	Test Cup Assembly for apm (with Handle)	1	
1000711	Thermofuse (Pack of 5, with Insulating tube)	1	
	Gas Hose $\varphi$ 15 x $\varphi$ 9.5, 1.5m	1	with Gas Hose Band
1007420	Electric Igniter, El-8 for abl-8, atg-8, apm-8, apm-100	1	
1006746	Windscreen for apm-100	1	
1006970	Screw for Windscreen for FP100	2	
1006922	AC Power Cord 100 V to 120 V 2 m	1	
1006930	AC Power Cord 220 V to 240 V 2 m	1	
	USB Flash Drive (Instruction & Maintenance Manual)	1	

## **Optional Accessories**

Code No.	Part Name	Qty	Remarks
1005593	Printer NEX-C200U01 AC 100 V to 120 V	1	
1005594	Printer NEX-C200U05 AC 220 V to 240 V	1	
1005861	Printer Roll Paper W58 $ imes$ $\varphi$ 50 for NEX-C200	1	

## **Suggested Spares for 2 years**

Code No.	Part Name	Qty	Remarks	
1007478	Test Cup Assembly for apm (with Handle)	1		
1000711	Thermofuse (Pack of 5, with Insulating tube)	1		
1007420	Electric Igniter, El-8 for abl-8, atg-8, apm-8, apm-100	5		
1001314	Temperature Sensor for apm-8, apm-100	1		
1001294	Flash Detector for APM-7, apm-8, apm-100	1		

Distributed by



**Tanaka Scientific Limited** 

7-10-3, Ayase , Adachi-ku, Tokyo,120-0005 JAPAN TEL:+81-3-3620-1711 E-mail:overseas-group@tanaka-sci.com FAX:+81-3-3620-1713 URL :www.tanaka-sci.com