

# Micro Carbon Residue Tester

Micro Carbon Residue Tester

acr-100m



## Connectable to Electronic Balance

Weight values can be imported and test results are calculated automatically.

## Equipped with Printing Function

Printing is available if a printer is connected.

## LCD Touch Screen

Equipped with a touch screen for simple, intuitive operation.

## Interlock Safety Function

Cannot be started unless the lid is in place.

## Cleaning the Oven Seal Lid Assembly

Possible to burn off the dirt from the lid assembly by leaving it in a coking oven.

## Test Results

Results from the last 6 test runs can be stored.



## Specifications



Test Methods	ASTM D4530, ISO 10370, JIS K2270-2
Measuring Range	Ambient to 500 °C
Heating / Cooling	Pipe Heater / Air Cooling
I/O Port	SD Card x 1, Ethernet x 1 (Printer / LIMS), RS-232C x1
Dimensions	W 390 mm x D 480 mm x H 480 mm
Weight	25 kg
Power Consumption	AC 100 V to 120 V 50 / 60 Hz 1500 VA AC 220 V to 240 V 50 / 60 Hz 1500 VA

## Tanaka Scientific Limited

7-10-3, Ayase, Adachi-ku, Tokyo, 120-0005 JAPAN  
TEL : +81-3-3620-1711 E-mail : overseas-group@tanaka-sci.com  
FAX : +81-3-3620-1713 URL : www.tanaka-sci.com